

## Carnét Pressure Regulator

### Safe, precise medical air therapy

This pressure regulator provides a safe and convenient way of using high pressure gas cylinders to connect to flowmeters, ventilators and other equipment required to run at low gas pressure.

#### Features

- State-of-the-art balanced piston design
- No membrane present
- Only approved polymers (halogen-free) in contact with the high pressure gas (as per ISO 15001)
- Long service interval

#### Benefits

- Improved stability of outlet pressure across inlet (cylinder) pressure range
- Reduces maintenance
- The use of halogen-free polymers significantly reduces the risk of patient exposure to toxic or explosive elements, should the regulator be exposed to high temperatures
- Low maintenance costs

#### Specifications

- Inlet pressure: up to 20,000 kPa (2,900 psi)
- Outlet pressure: 400kPa (58 psi)
- Inlet fitting: pin index yoke
- Outlet fitting: AS sleeve index system or BS schrader

#### Manufacturer

BPR Medical

#### Warranty Period

1 year

#### Service Interval

4 years



Outlet: AS



Outlet: BS

### Quality Assurance

This device conforms to the following standards:

- ISO 10524.1: Pressure regulators for use with medical gases
- ISO 15001: Anaesthetic and respiratory equipment
- AS 3840.1: Pressure regulators for use with medical gases
- ISO 15024.1: Fibre-reinforced plastic composites
- AS2896
- BS5682

### Pressure Regulator

Part Number **AS**  
**BS**

819-0077 AUST. AND N.Z.

819-0082 N.Z. ONLY

To find out more or place an order contact our customer service team:

Australia

**1800 050 999**

or visit [www.boc-healthcare.com.au](http://www.boc-healthcare.com.au)

New Zealand

**0800 656 334**

or visit [www.boc-healthcare.co.nz](http://www.boc-healthcare.co.nz)

**Australia: BOC Limited** ABN 95 000 029 729, BOC Healthcare 10 Julius Ave, North Ryde NSW 2113

T: 1300 363 109 F: 1300 363 438 E: [hospital.care@voc.com](mailto:hospital.care@voc.com)

**New Zealand: BOC Limited** WN007748, BOC Healthcare NZ, 988 Great South Rd, Penrose Auckland,

T: 0800 656 334 F: 0800 275 275 E: [HealthcareNZ@voc.com](mailto:HealthcareNZ@voc.com)

Details given in this document are believed to be correct at the time of printing. While proper care has been taken in the preparation, no liability for injury or damage resulting from its use can be accepted.

BOC is a trading name of BOC Limited, a member of The Linde Group. © BOC Limited 2013.

Reproduction without permission is strictly prohibited.